





Endobrow Tungsten Microdissection Needles



ASSI.37123
 7" long
 (18 cm)

ASSI.36123
 6" long
 (15 cm)



- Ultra sharp Tungsten needle tips
- Cut and dissect tissue cleanly, with minimal blood loss or tissue trauma
- Focused Monopolar energy at the needle tip
- Non-sterile packed for single use
- Insulated shaft 3/32" (0.24cm) sleeve 5/8" (1.59cm)

Endobrow
ASSI.36123 6" (15 cm)
ASSI.37123 7" (18 cm)