

# ASSI® Rake Retractors & THIMBLE HOOKS

RAKE RETRACTORS • THIMBLE HOOKS

ASSI.AG32526



ASSI.AG33526



ASSI.AG32226



ASSI.AG32126



ASSI.AG31426



ASSI.AG430526



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®

*For diamond perfect performance®*

accurate surgical & scientific instruments corporation

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

[www.accuratesurgical.com](http://www.accuratesurgical.com)

SEE COMPLETE LINE  
OF PRODUCTS ON  
THE OTHER SIDE



## WHAT'S NEW

### ASSI.AG383326


Gorney Maxwell Rake Retractor with spring,  
3 1/2" (7.5cm) 38mm wide



## ABDOMINAL RETRACTOR

### ASSI.AG32526

Abdominal Rake Retractor, 8 1/4" (21cm),  
double ended, 6 sharp prongs,  
(15cm) wide, s.s.



## CORONAL RETRACTORS


### ASSI.AG68826

Coronal Lift Retractor,  
8 1/4" (21cm), 6 sharp prongs,  
6" (15cm) wide, s.s.




### ASSI.AG32126

Coronal Brow Lift Retractor, 7" (18cm), double ended,  
7 sharp prongs, 2 3/4" (7cm) wide, s.s.



### ASSI.AG32226R

Coronal Brow Lift Retractor, 7" (18cm), double ended,  
7 sharp prongs, reverse, 2 3/4" (7cm) wide, s.s.



## THIMBLE HOOKS

### ASSI.AG430426

Millard Type Thimble Hook, single hook



### ASSI.AG430526

Millard Type Thimble Hook, double hook



## RAKE RETRACTORS

### ASSI.AG31426

Freeman Rake Retractor, 7" (18cm), double ended,  
4 sharp prongs, 2 recessed, 1 1/2" (4cm) wide, s.s.



### ASSI.AG38326

Freeman Rake Retractor, 7" (18cm), double ended,  
4 sharp prongs, all same length, 1 1/2" (4cm) wide, s.s.



### ASSI.AG31926

Freeman Rake Retractor, 7" (18cm), double ended,  
4 sharp prongs, 2 recessed, 7/8" (2.2cm) wide, s.s.



### ASSI.AG31826

Freeman Rake Retractor, 7" (18cm), double ended,  
4 sharp prongs, all same length, 7/8" (2.2cm) wide, s.s.




### ASSI.AG31726

ASSI® Mini Rake Retractor, 6 1/8" (15.5cm),  
double ended, 4 semi sharp prongs, 3/4" (2cm) wide, s.s.



### ASSI.AG13626

Freeman Adamson Face Lift Retractor, 7" (18cm),  
4 sharp prongs, angled, 1 1/2" (4cm) wide, s.s.



### ASSI.AG312126

Freeman Rake Retractor, 2" (5cm), with thumb  
ring, 4 sharp prongs, 3/4" (2cm) wide, s.s.




### ASSI.AG419126

Same as above with Reversed Prongs.

### ASSI.AG32226

Freeman Rake Retractor, 2" (5cm), short  
handle with thumb ring, 4 sharp prongs,  
1 1/2" (4cm) wide, s.s.



### ASSI.AG519126

Same as above with Reversed Prongs.

### ASSI.AG31526

Freeman Rake Retractor, 4" (10cm), short  
handle with thumb ring, 4 sharp prongs,  
1 1/2" (4cm) wide, s.s.



### ASSI.AG33526

Rake Retractor/Skin Hook, (11.5cm), 3 sharp prongs, s.s.



**assi**®

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®

For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

[www.accuratesurgical.com](http://www.accuratesurgical.com)