

# EPSTEIN ABDOMINOPLASTY RETRACTOR

As the trend in today's abdominoplasties is towards less upper lateral undermining, the upper tunnel has become narrower. During the dissection of this tunnel, it is difficult to maintain countertraction with available instrumentation. The Epstein Abdominoplasty Retractor was designed to assist in the performance of the dissection of the upper abdominal tunnel. The ergonomic handle is easily held by the surgical assistant. It is available in several blade lengths so as to best fit the anatomy of the patient: whether the tunnel is long or short, there is a retractor to provide the best mechanical advantage in yielding exposure and reducing fatigue. The widened, curved working end spreads the tissues of the upper skin flap apart as the teeth gently hold them in place without slippage. The leading edge of the dissection is easily seen and maintained, facilitating effortless cautery elevation of the skin flap from the muscle fascia. The retractor is also extremely helpful in elevating the abdominal skin flap over the narrow tunnel so that the underlying muscle fascia can be plicated.

- Designed to assist in the performance of the dissection of the upper abdominal tunnel
- The wide curved working end spreads and holds the tissue of the upper skin flap apart
- Facilitates cautery elevation of the skin flap from the muscle fascia
- Extremely helpful in elevating the skin flap over the narrow tunnel for muscle fascia plication

#### **ASSI.ABR31726**

#### **Epstein Abdominoplasty Retractor**

80x27mm blade with 60mm wide working end with teeth and ergonomic handle

#### **ASSI.ABR31826**

#### **Epstein Abdominoplasty Retractor**

110x27mm blade with 60mm wide working end with teeth and ergonomic handle

#### **ASSI.ABR32026**

#### **Epstein Abdominoplasty Retractor**

140x27mm blade with 60mm wide working end with teeth and ergonomic handle

#### **ASSI.ABR31926**

#### **Epstein Abdominoplasty Retractor**

160x27mm blade with 60mm wide working end with teeth and ergonomic handle

#### **ASSI.ABR32226**

#### **Epstein Abdominoplasty Retractor**

240x27mm blade with 60mm wide working end with teeth and ergonomic handle

**assi**®   
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
For diamond perfect performance®

**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

[www.accuratesurgical.com](http://www.accuratesurgical.com)

Designed By:  
Mark D. Epstein,  
M.D., F.A.C.S.,  
Stony Brook, NY