


# A Cool Dry Alternative to Messy Wet Gauze Pads



## Gel Masks and Eye Patches from ASSI.®

- Sterile blue liquid that retains cold or heat for up to one hour
- May be kept refrigerated
- Can be warmed in a pan of hot water to desired temperature
- Easy to clean and reusable
- Elastic nylon band to hold the mask in place
- Designed so the patches rest comfortably in the eye socket



**ASSI.A42343**  
Gel Eye Mask, 10 per box



**ASSI.A43343**  
Gel Eye Patches, 10 pair per box



*For diamond perfect performance®*

**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

e-mail: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

internet: [www.accuratesurgical.com](http://www.accuratesurgical.com)