

# ASSI® Double Handle Breast Retractors



The ASSI® Double Handle Breast Retractors are designed so both the surgeon and the assistant can grasp the retractor during Breast Augmentation Pocket Dissection, thus alleviating surgeon fatigue.

They are available in different blade lengths and widths, with and without teeth as well as with the suction tube in the standard position or reversed 180°.



ASSI.ABR55726



ASSI.ABR55926  
w/oval shape fiber optic



ASSI.ABR54826  
w/teeth

Designed by: Mark A. Schusterman, M.D. FACS; Houston, TX

## ASSI® Breast Retractor With Double Handle

**ASSI.ABR54826**  
150x25mm, w/teeth,  
w/fiber optic, suction & adapter  
**w/reverse suction tube**

**ASSI.ABR54926**  
150x25mm, w/teeth,  
w/fiber optic & adapter

**ASSI.ABR55026**  
80x25mm, w/teeth,  
w/fiber optic, suction & adapter  
**w/reverse suction tube**

**ASSI.ABR55126**  
80x25mm, w/teeth,  
w/fiber optic, & adapter

**ASSI.ABR55226**  
180x25mm,  
w/fiber optic, suction & adapter  
**w/reverse suction tube**

**ASSI.ABR55426**  
150x16mm,  
w/fiber optic, suction & adapter

**ASSI.ABR55526**  
180x16mm,  
w/fiber optic, suction & adapter

**ASSI.ABR55926**  
150x30mm,  
w/oval shape fiber optic, suction &  
teeth

**ASSI.ABR55726**  
150x25mm,  
w/fiber optic, suction & adapter

**ASSI.ABR56026**  
80x16mm,  
w/fiber optic, suction & adapter  
**w/reverse suction tube**

**ASSI.ABR56426**  
150x16mm,  
w/fiber optic, suction & adapter  
**w/reverse suction tube**

**ASSI.ABR56526**  
180x16mm,  
w/fiber optic, suction & adapter  
**w/reverse suction tube**

**ASSI.ABR56726**  
180x16mm,  
w/fiber optic, suction & adapter  
**w/teeth**

**ASSI.ABR56826**  
180x30mm,  
w/fiber optic, suction & adapter



**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

[www.accuratesurgical.com](http://www.accuratesurgical.com)