

assi[®]

NOBLE[™]

REAL NON-STICK BIPOLAR FORCEPS & MICRO MONOPOLAR FORCEPS



Now
Available With
Irrigation!

- The key to the Noble[™] Metal's Non-Stick property is that the energy is moved to and from the distal tips and contact surfaces many times faster than conventional Stainless Steel.
- The rapid energy movement eliminates "hot spots" and prevents tissue buildup on the forceps.
- The rapid energy movement results in lower temperature settings: No "charring" and minimizes "scarring".
- The Noble[™] Metal Non-Stick jaw is permanent. It will not be compromised by frequent sterilizations. The jaws are made from the Noble Metal, not just a coating or highly polished surface.
- The Noble[™] Metal Non-Stick Forceps are extremely economical: They save the surgeon valuable time and frustration as the tips do not continuously have to be removed from the operative site to be cleaned.

assi[®] 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS[®]

For diamond perfect performance[®]

accurate surgical & scientific instruments corporation


300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

Not all ASSI products shown in our literature or on our website are available for sale in Canada.


NON-STICK BIPOLAR FORCEPS



Size	2-Pin Round Plug Product Number	Flat Plug (not shown) Product Number	Tip Diameter	Tip Style
Jewelers				
11cm	ASSI.BPNS11223	ASSI.BPNS11223FP	0.25mm	Straight, Pointed Tip
	ASSI.BPNS11224	ASSI.BPNS11224FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS11226	ASSI.BPNS11226FP	1.0mm	Straight, Blunt Tip
11cm	ASSI.BPNS11233	ASSI.BPNS11233FP	0.25mm	Angle, Pointed Tip
	ASSI.BPNS22467	ASSI.BPNS22467FP	0.6mm	Angle, Blunt Tip
	ASSI.BPNS11236	ASSI.BPNS11236FP	1.0mm	Angle, Blunt Tip
Adson:				
12cm	ASSI.BPNS12723		0.3mm	Straight
	ASSI.BPNS12729	ASSI.BPNS12729FP	0.6mm	Straight
	ASSI.BPNS12724	ASSI.BPNS12724FP	1.0mm	Straight
	ASSI.BPNS12727		1.2mm	Straight
	ASSI.BPNS12726	ASSI.BPNS12726FP	1.5mm	Straight
Bayonet:				
18cm	ASSI.BPNS18423	ASSI.BPNS18423FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS18429	ASSI.BPNS18429FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS18424	ASSI.BPNS18424FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS18425	ASSI.BPNS18425FP	2.0mm	Straight, Blunt Tip
20cm	ASSI.BPNS20423	ASSI.BPNS20423FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS20429	ASSI.BPNS20429FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS20424	ASSI.BPNS20424FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS20426		1.2mm	Straight, Blunt Tip
	ASSI.BPNS20425	ASSI.BPNS20425FP	2.0mm	Straight, Blunt Tip
	ASSI.BPNS20454	ASSI.BPNS20454FP	1.0mm	Angle Up, Blunt Tip
22cm	ASSI.BPNS22423	ASSI.BPNS22423FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS22429	ASSI.BPNS22429FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS22424	ASSI.BPNS22424FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS22425	ASSI.BPNS22425FP	2.0mm	Straight, Blunt Tip
	ASSI.BPNS22454	ASSI.BPNS22454FP	1.0mm	Angle Up, Blunt Tip
	ASSI.BPNS22464	ASSI.BPNS22464FP	1.0mm	Angle Down, Blunt Tip
24cm	ASSI.BPNS24423	ASSI.BPNS24423FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS24429	ASSI.BPNS24429FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS24424	ASSI.BPNS24424FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS24426		1.2mm	Straight, Blunt Tip
	ASSI.BPNS24425	ASSI.BPNS24425FP	2.0mm	Straight, Blunt Tip




accurate surgical & scientific instruments corporation
 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
 800.645.3569 • Info: assi@accuratesurgical.com
 Orders: orders@accuratesurgical.com
www.accuratesurgical.com



Size	2-Pin Round Plug Product Number	Flat Plug (not shown) Product Number	Tip Diameter	Tip Style
Straight				
15cm	ASSI.BPNS15323	ASSI.BPNS15323FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS15329	ASSI.BPNS15329FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS15324	ASSI.BPNS15324FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS15325	ASSI.BPNS15325FP	2.0mm	Straight, Blunt Tip
15cm	ASSI.BPNS15333	ASSI.BPNS15333FP	0.3mm	Angled, Pointed Tip
	ASSI.BPNS15339	ASSI.BPNS15339FP	0.6mm	Angled, Blunt Tip
	ASSI.BPNS15334	ASSI.BPNS15334FP	1.0mm	Angled, Blunt Tip
	ASSI.BPNS15335	ASSI.BPNS15335FP	2.0mm	Angled, Blunt Tip
18cm	ASSI.BPNS18323	ASSI.BPNS18323FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS18329	ASSI.BPNS18329FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS18324	ASSI.BPNS18324FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS18325	ASSI.BPNS18325FP	2.0mm	Straight, Blunt Tip
18cm	ASSI.BPNS18333	ASSI.BPNS18333FP	0.3mm	Angled, Pointed Tip
	ASSI.BPNS18339	ASSI.BPNS18339FP	0.6mm	Angled, Blunt Tip
	ASSI.BPNS18334	ASSI.BPNS18334FP	1.0mm	Angled, Blunt Tip
	ASSI.BPNS18335	ASSI.BPNS18335FP	2.0mm	Angled, Blunt Tip
	ASSI.BPNS18395	ASSI.BPNS18395FP	2.0mm	Straight, Blunt Serrated Tip
20cm	ASSI.BPNS20323	ASSI.BPNS20323FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS20329	ASSI.BPNS20329FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS20324	ASSI.BPNS20324FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS20325	ASSI.BPNS20325FP	2.0mm	Straight, Blunt Tip
22cm	ASSI.BPNS22323	ASSI.BPNS22323FP	0.3mm	Straight, Pointed Tip
	ASSI.BPNS22329	ASSI.BPNS22329FP	0.6mm	Straight, Blunt Tip
	ASSI.BPNS22324	ASSI.BPNS22324FP	1.0mm	Straight, Blunt Tip
	ASSI.BPNS22325	ASSI.BPNS22325FP	2.0mm	Straight, Blunt Tip
24cm	ASSI.BPNS24325		2.0mm	Straight, Blunt Tip
26cm	ASSI.BPNS26325		2.0mm	Straight, Blunt Tip
28cm	ASSI.BPNS28325		2.0mm	Straight, Blunt Tip

Bayonet, With Irrigation				
20cm	ASSI.BPNS20423IR		0.3mm	Straight, Pointed Tip
	ASSI.BPNS20424IR		1.0mm	Straight, Blunt Tip
	ASSI.BPNS20425IR		2.0mm	Straight, Blunt Tip
	ASSI.BPNS20429IR		0.6mm	Straight, Blunt Tip




NON-STICK BIPOLAR FORCEPS




accurate surgical & scientific instruments corporation
 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
 800.645.3569 • Info: assi@accuratesurgical.com
 Orders: orders@accuratesurgical.com
www.accuratesurgical.com

RACKS



ASSI.72259
Rack for Bayonet Bipolar Forceps
2 finger & 1 clamping bar
Large, 25cm x 25cm



ASSI.72257
Rack for Bipolar Forceps/Bipolar Scissors
2 finger & 1 clamping bar
Small, 23cm x 22cm

POLISHING CLOTH



Ultra Soft Polishing Cloth Two – Cloth System in 1 Cloth

ASSI.NBPPC46
Polishing Cloth


- The ASSI® Noble™ Bipolar Forceps Polishing Cloth System™ may be used to remove oxidation from the Noble™ Metal Tip.
- After multiple reuse and sterilization cycles, a layer of oxidation may form on the surface of the Noble™ Bipolar Forceps Tips. Oxidation on the tip is harmless and does not interfere with the non-stick performance.
- Once you have decontaminated and cleaned the forceps use the Inner White Cloth to remove the oxidation.
- Use the Outer Blue Cloth to gently polish the tips before you autoclave.
- Repeat this process when and if oxidation returns.

ASSI.NBPPC46 Noble Bipolar Forceps Polishing Cloth: 3 each 4 x 6" cloths per package

NOBLE™ MICRO MONOPOLAR NON-STICK FORCEPS

Insulated, hand-activated with **Non-Stick** straight, smooth jaws.
Eliminates “hot pots” and prevents tissue build up on the forceps


- The key to the Noble™ Metal's Non-Stick property is the Rapid energy movement resulting in Lower Temperature Settings: No Charring and Minimizes Scarring.
- The Noble™ Metal Non-Stick jaw is Permanent made from the Noble Metal. Not just a coating or highly polished surface.
- The tips do not continuously have to be removed from the operative site to be cleaned.
- ASSI® Noble™ Micro Monopolar Forceps are specifically designed to reach the entire length of the pocket.
- Ideal for Inframammary and Initial, Upper Axillary Pocket Dissections.
- The forceps precisely combines cut and coagulation current for best possible dissection averting current spread, tissue damage and bleeding.



601240MMPNS
Non-Stick Bayonet Forceps
24cm, Insul. Hnd. Actv.Str. Jaw
0.7mm Tip
Limited Reuseable

601900MMPNS
Non-Stick Bayonet Forceps
19cm, Insul. Hnd. Actv. Str. Jaw
0.7mm Tip
Limited Reuseable

601500MMPNS
Non-Stick Bayonet Forceps
15cm, Insul. Hnd. Actv. Str. Jaw
0.7mm Tip
Limited Reuseable



Special Non-Stick
Silver Jaws




602400MMPNS
Non-Stick Straight Forceps
24cm, Insul. Hnd. Actv.
Str. Jaw 0.7mm Tip
Limited Reuseable

601270MMPNS
Non-Stick Straight Forceps
26cm, Insul. Hnd. Actv.
Str. Jaw 0.7mm Tip
Limited Reuseable


assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation
516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
800.645.3569 • Info: assi@accuratesurgical.com
Orders: orders@accuratesurgical.com
www.accuratesurgical.com


COAGULATING CORDS




Rubber Bipolar Cord
Female end will fit standard two-pin forceps connector.
228000C



Rubber Bipolar Cord With 45° Angle Forceps Plug
Female end will fit standard two-pin forceps connector.
228001C



Silicone Bipolar Cord With 45° Angle Forceps Plug
Female end will fit standard two-pin forceps connector.
229001C



Disposable Sterile Bipolar Cord
(One time use only)
Female end will fit standard two-pin forceps connector.
227000C



ASSI.1766 Cable Silicone For Bipolar Forceps Non-Stick
10 ft., 2 pin plug w/ Valley Lab 2 prong plug

ASSI.1788 Cable Silicone For Bipolar Forceps Non-Stick
17 ft., 2 pin plug w/ Valley Lab 2 prong plug



FOR HANDSWITCH FORCEPS
SILICONE BIPOLAR HANDSWITCH CORD
Female end will fit standard three-pin forceps connector.
230001C




For Noble Non-Stick Bipolar Forceps With Flat Plug
Silicone Bipolar Cord F/ Flat Plug Bipolar Forceps w/ Banana Plugs, 10 ft.
221010C




ASSI.2766 Cable Silicone For Bipolar Forceps Non-Stick
10 ft., 2 pin angled plug w/ Valley Lab 2 prong plug

ASSI.2788 Cable Silicone For Bipolar Forceps Non-Stick
17 ft., 2 pin angled plug w/ Valley Lab 2 prong plug



ASSI.1566 Silicone Cord For Bipolar Forceps Non-Stick
17ft, Square flat plug w/ Valley Lab 2 prong plug

ASSI.1577 Silicone Cord For Bipolar Forceps Non-Stick
10ft, Square flat plug w/ Valley Lab 2 prong plug



ASSI.3422 Silicone Cord For Bipolar Forceps Non-Stick
10ft, Square flat plug w/ Valley Lab 3 prong plug