

# ASSI® Forehead and Face Lift Instruments are like Diamonds...



*Crafted for Perfection.*

*Created for Performance.*

*Cut with Precision...*

*the way you do.*

**Available**  
with Ergonomic  
Handles

ASSI.AEP222326E

ASSI.AEP14826

**assi**® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
*For diamond perfect performance*®


accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.AEP715236




**ASSI.AEP715236**  
Scalp Elevator, 24cm,  
half curve, **12.5mm wide blade**  
with beveled sharp edges

ASSI.AEP715236E



**ASSI.AEP715236E**  
with ergonomic handles

ASSI.AEP715526



**ASSI.AEP715526**  
Scalp Elevator 24cm,  
half curve, **12.5mm wide blade**  
with beveled sharp edges

ASSI.AEP715526E



**ASSI.AEP715526E**  
with ergonomic handles

ASSI.AEP715726




**ASSI.AEP715726**  
Scalp Elevator dissector, 24cm,  
straight

ASSI.AEP715726E



**ASSI.AEP715726E**  
with ergonomic handles

ASSI.AEP112726




**ASSI.AEP112726**  
Scalp Elevator, 24cm  
"S"-shaped , shaft bent 10° at  
handle, quarter curve, **12.5mm**  
**wide blade** with beveled sharp  
edges

ASSI.AEP112726E



**ASSI.AEP112726E**  
with ergonomic handles

ASSI.AEP715926



**ASSI.AEP715926**  
Scalp Elevator, 24cm,  
strong half curve,  
**15mm wide blade**  
with beveled sharp edges

ASSI.AEP715926E



**ASSI.AEP715926E**  
with ergonomic handles




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378


Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.AEP723926



ASSI.AEP723926E



**ASSI.AEP723926**  
Scalp elevator, 43cm;  
32cm insulated shaft,  
strong half curve,  
15mm wide blade with  
beveled sharp edges

**ASSI.AEP723926E**  
with ergonomic handles

ASSI.AEP715826



ASSI.AEP715826E



**ASSI.AEP715826**  
Scalp Elevator Dissector, 24cm,  
slightly curved

**ASSI.AEP715826E**  
with ergonomic handles

ASSI.AEP715626




ASSI.AEP715626E



**ASSI.AEP715626**  
Scalp Elevator Dissector, 24cm,  
half curve

**ASSI.AEP715626E**  
with ergonomic handles

ASSI.AEP725926




ASSI.AEP725926E




**ASSI.AEP725926**  
Scalp Elevator Dissector, 24cm,  
curved, 7mm wide

**ASSI.AEP725926E**  
with ergonomic handles

ASSI.AEP714826



Not available with Ergonomic Handle



**ASSI.AEP714826**  
Scalp Elevator, 26cm,  
slightly curved, with suction

ASSI.AEP727626



Not available with Ergonomic Handle



**ASSI.AEP727626**  
Scalp Elevator Dissector, 24cm,  
curved, 7mm wide blade  
with suction



accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.AEP212326



**ASSI.AEP212326**  
Frontotemporal Dissector, 22.5cm,  
curved, with beveled edges

ASSI.AEP212326E



**ASSI.AEP212326E**  
with ergonomic handles

ASSI.AEP14326




**ASSI.AEP14326**  
Frontotemporal Dissector, 23cm,  
straight, **10mm wide blade**  
with beveled sharp edges

ASSI.AEP14326E



**ASSI.AEP14326E**  
with ergonomic handles

ASSI.AEP232326




**ASSI.AEP232326**  
Temporal T-Dissector, 23.5cm,  
straight

ASSI.AEP232326E



**ASSI.AEP232326E**  
with ergonomic handles

ASSI.AEP222326




**ASSI.AEP222326**  
Parietal Periosteal elevator, 25cm,  
curved

ASSI.AEP222326E



**ASSI.AEP222326E**  
with ergonomic handles

ASSI.AEP14826



**ASSI.AEP14826**  
Periosteal Spreader, 24cm,  
**8mm wide "S"-shaped blade**

ASSI.AEP14826E



**ASSI.AEP14826E**  
with ergonomic handles




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.AEP38326




**ASSI.AEP38326**  
Periosteal Spreader, 24cm;  
10mm wide "S"-shaped blade

ASSI.AEP38326E



**ASSI.AEP38326E**  
with ergonomic handles

ASSI.AEP242326




**ASSI.AEP242326**  
Frontoglabellar Dissector,  
26.5cm, "S"-shaped curved  
with beveled edges

ASSI.AEP242326E



**ASSI.AEP242326E**  
with ergonomic handles

ASSI.AEP252326




**ASSI.AEP252326**  
Arcus Marginalis Dissector, 24cm,  
"S"-shaped curved

ASSI.AEP252326E



**ASSI.AEP252326E**  
with ergonomic handles

ASSI.AEP272326




**ASSI.AEP272326**  
Flap Dissector, 23.5cm,  
straight with beveled edges

ASSI.AEP272326E



**ASSI.AEP272326E**  
with ergonomic handles

ASSI.AEP282326



**ASSI.AEP282326**  
Zygomatic Arch Elevator, 20.5cm,  
straight with beveled edges

ASSI.AEP282326E



**ASSI.AEP282326E**  
with ergonomic handles




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.AEP292326




**ASSI.AEP292326**  
Midface Facial Dissector, 22.5cm,  
half curve

ASSI.AEP292326E



**ASSI.AEP292326E**  
with ergonomic handles

ASSI.AEP716226




**ASSI.AEP716226**  
Orbital Rim Dissector, 24cm,  
half curve

ASSI.AEP716226E



**ASSI.AEP716226E**  
with ergonomic handles

ASSI.AEP712426




**ASSI.AEP712426**  
Transoral Dissector, 23.5cm,  
curved

ASSI.AEP712426E



**ASSI.AEP712426E**  
with ergonomic handles

ASSI.AEP712826



**ASSI.AEP712826**  
Corrugator and Procerus Muscle  
Dissector, 19.5c

Not available with Ergonomic Handle




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada


ASSI.AEP712526



Not available with Ergonomic Handle

**ASSI.AEP712526**  
Nerve Protector, 17cm **curved left**


ASSI.AEP712626



Not available with Ergonomic Handle

**ASSI.AEP712626**  
Nerve Protector, 17cm **curved right**

ASSI.AEP713626



**ASSI.AEP713626**  
Nerve Cutting Knife, 17cm  
3mm blade

ASSI.AEP724126 LEFT ANGLE HOOK



ASSI.AEP724126E



ASSI.AEP724236 RIGHT ANGLE HOOK



ASSI.AEP724236E



**ASSI.AEP724126**  
Nerve Hook, 24cm, half curve, finely tapered, **left angle hook**

**ASSI.AEP724126E**  
with ergonomic handles

**ASSI.AEP724236**  
Nerve Hook, 24cm, half curve, finely tapered, **right angle hook**

**ASSI.AEP724236E**  
with ergonomic handles

ASSI.AEP716126 RIGID SHAFT



ASSI.AEP716926 MALLEABLE SHAFT



Not available with Ergonomic Handle

**ASSI.AEP716126**  
Nerve Hook, **rigid shaft**, 24cm, straight

**ASSI.AEP716926**  
Nerve Hook, **malleable shaft**, 24cm, straight




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.AEP716026




**ASSI.AEP716026**  
Nerve Dissector, 24cm,  
quarter curve

ASSI.AEP716026E




**ASSI.AEP716026E**  
with ergonomic handles

ASSI.AEP714236




**ASSI.AEP714236**  
Nerve Dissector, 24cm,  
slightly curved, spatulated tip

ASSI.AEP714236E




**ASSI.AEP714236E**  
with ergonomic handles

ASSI.AEP727236



**ASSI.AEP727236**  
Knife Handle for plastic surgery,  
20.5cm, slight curve

Not available with Ergonomic Handle



**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI.SU263326



Sheaths with **rotating** and **fixed** stop-cocks, 17.5cm, **4mm dia.**; with non-glare satin finish.

**ASSI.SU263326**

**1 rotating** stop-cock

**ASSI.SU52426**

**1 fixed** stop-cock (not shown)

ASSI.SU46326



**ASSI.SU46326**

**2 rotating** stop-cocks

**ASSI.SU36326**

**2 fixed** stop-cocks (not shown)

ASSI.SU374526



Sheaths with **rotating** stop-cock, 17.5cm, **5mm dia.**; with non-glare satin finish.

**ASSI.SU374526**

**1 rotating** stop-cock

ASSI.SU62426



Sheath 100 mm working length, **4mm dia**, **rotating** stop cock with angled spatulated tip.

**ASSI.SU62426**

**1 rotating** stop-cock

**ASSI.SU62726**

**2 rotating** stop-cocks (not shown)

Sheath 100 mm working length, **5mm dia**, **rotating** stop cock with angled spatulated tip.

**ASSI.SU62626**

**1 rotating** stop-cock (not shown)

**ASSI.SU62826**

**2 rotating** stop-cocks (not shown)

ASSI.SU63026



Sheath 100 mm working length, **4mm dia**, **rotating** stop cock with **wide** angled spatulated tip.

**ASSI.SU63026**

**1 rotating** stop-cock

Sheath 100 mm working length, **5mm dia**, **rotating** stop cock with **wide** angled spatulated tip.

**ASSI.SU63226**

**1 rotating** stop-cock (not shown)




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada



ASSI.AEP216326

**SCISSORS, 12.5cm**

ASSI.AEP216026 curved down	ASSI.AEP217126 23cm, left curved, 2.7mm dia.	ASSI.AEP217226 23cm, right curved, 2.7mm dia.
ASSI.AEP216126 left curved	ASSI.AEP216226 right curved	ASSI.AEP216326 straight



ASSI.AEP383126




ASSI.AEP383226



ASSI.AEP383326

**SCISSORS, 15cm, Ergonomic Handle (Ramirez Type)**

ASSI.AEP383126 straight malleable shaft, blade curved left	ASSI.AEP383226 curved left malleable shaft, straight blade	ASSI.AEP383326 curved right malleable shaft, straight blade
--	--	---



ASSI.AEP216626

**FORCEPS, 12.5cm, Grasping**

ASSI.AEP216626 left curved	ASSI.AEP216726 right curved	ASSI.AEP216926 straight
-------------------------------	--------------------------------	----------------------------



ASSI.AEP383426



ASSI.AEP383526

**FORCEPS, 15cm, Ergonomic Handle (Ramirez Type), straight malleable shaft**

ASSI.AEP383426 Biopsy	ASSI.AEP383526 Grasping
--------------------------	----------------------------



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®


For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada



ASSI.AEP42326




**FORCEPS**

12.5cm, Grasping, Insulated Shaft, Double Action Jaw, 2.7mm dia.

ASSI.AEP42326  
30° left curved  
ASSI.AEP52326  
30° right curved

ASSI.AEP75026  
serrated straight  
ASSI.AEP33826  
30° angled downward

ASSI.AEP216826



**PUNCH, 12.5cm**

ASSI.AEP216426  
right curved

ASSI.AEP216526  
left curved

ASSI.AEP216826  
straight

ASSI.AEP383626



**NEEDLE HOLDER TC, 15cm, Ergonomic Handle (Ramirez Type)**

ASSI.AEP383626 straight malleable shaft

**assi**<sup>®</sup>



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS<sup>®</sup>

*For diamond perfect performance*<sup>®</sup>

**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada



**ASSI.AEP73326I**  
Micro Maryland Dissecting Forceps/Grasper  
cvd shaft, rotating, insulated,  
Insert Only



**ASSI.AEP72726I**  
Micro Grasping Forceps  
cvd shaft, rotating, insulated,  
Insert Only



**ASSI.AEP73126I**  
Micro Grasping Forceps  
Alligator Jaw, cvd shaft, rotating, insulated,  
Insert Only



**ASSI.AEP72326I**  
Micro Hook Scissors  
cvd shaft, rotating, insulated,  
Insert Only




**ASSI.AEP73726I**  
Micro Metzenbaum Scissors  
cvd blade, cvd shaft, rotating, insulated,  
Insert Only



**ASSI.AEP73826I**  
Micro Metzenbaum Scissors  
str. blade, cvd shaft, rotating, insulated,  
Insert Only



**ASSI.AEP73426I**  
Micro Needle Holder Insert,  
working length 15 cm. 3.0 mm shaft dia.  
single action rotating, insulated,  
Insert Only




**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

**ASSI.AEP87326**  
Carbon Fiber Handle  
for Modular Instrument Inserts,  
rotatable,  
w/HF-connection



\*Handle for use with inserts  
on previous page

**ASSI.AEP87326R**  
Carbon Fiber Handle  
for Modular Instrument Inserts,  
rotating,  
w/ratchet w/HF-connection



\*Handle for use with inserts  
on previous page



**ASSI.AEP87426C**  
Cable for Carbon Fiber Handle  
3M long, 10ft. long for instruments  
w/4mm Valley Lab plug



**ASSI.AEP87726C**  
Cable for Carbon Fiber Handle  
3M long, 10ft. long  
w/4mm Banana plug

**ASSI.AEP87526C**  
Cable for Carbon Fiber Handle  
5M long, 15ft. long, for instruments  
w/4mm Valley Lab plug

**ASSI.AEP87926C**  
Cable for Carbon Fiber Handle  
5M long, 15ft. long  
w/4mm Banana plug



**ASSI.AEP112226**  
Cleaning Adapter with Luer-Lock  
for cleaning Endo Forehead &  
Facelift Instruments



**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

# Accessories



## Surgical Instrument Color Coding Tapes

ASSI® Color Coding Tapes make it easy to create your own labeling system for identifying and sorting instruments by type, specialty, even ownership.

These durable tapes come in 9 high-visibility colors. They are easy to apply. They have been clinically tested and found safe for use in conventional, flash and gas sterilization units.

ASSI Color Coding Tapes are 1/4" wide and packaged on spools in 300" lengths. Enough tape to color code approximately 230 instruments.

<b>ASSI.823</b>	Green	
<b>ASSI.825</b>	Yellow	
<b>ASSI.827</b>	Blue	
<b>ASSI.829</b>	Red	
<b>ASSI.831</b>	Brown	
<b>ASSI.833</b>	White	
<b>ASSI.835</b>	Orange	
<b>ASSI.837</b>	Black	
<b>ASSI.839</b>	Purple	

## Tip Covers– to prevent damage to instrument tips and edges

*Now available vented or non-vented*



- 11 Different sizes, each one color coded
- Autoclavable and steam-permeable
- Vented covers aid in the sterilizing and drying of condensation buildup on tips
- Stretchable — conforms to various instrument tip configurations
- Disposable
- Last for months
- Can be used on scissors, forceps, chisels, microsurgical and laparoscopic instruments, hemostats, curettes, rasps, knives, osteotomes, gouges, hooks, calipers, etc.
- Can be cut to any length
- Non-sterile

### ATC200 Mixed Assortment

(100 pieces, equal number of all sizes except ATC111-L)

### ATC200-V Mixed Assortment-Vented

(100 pieces, equal number of all sizes except ATC117 and ATC111-L)

### ATC400 Microsurgical Assortment

(100 pieces - 25 each: ATC100, 110, 111 and 112)

### ATC400-V Microsurgical Assortment-Vented

(100 pieces - 25 each: ATC100-V, 110-V, 111-V and 112-V)

### ATCBG Select Microsurgical Assortment

(100 pieces - 50 each: ATC100 and ATC110)

### ATCBG-V Select Microsurgical Assortment-Vented

(100 pieces - 50 each: ATC100-V and ATC110-V)

**Individual Package** (100 pieces, order by product number & color)  
**Non-vented**  
**Vented**

## VENTED AND NON-VENTED SIZES AND COLORS


Tip covers shown actual sizes

### NON-VENTED

ATC000    ATC100    ATC110    ATC111    ATC112    ATC113    ATC114    ATC115    ATC116    ATC117    ATC111-L

### VENTED

ATC000-V    ATC100-V    ATC110-V    ATC111-V    ATC112-V    ATC113-V    ATC114-V    ATC115-V    ATC116-V



Sold separately  
100 piece pkgs.

Not available in vented

# **assi**<sup>®</sup> NEW G7 CATALOG

*Be sure to ask for  
your copy of our new  
"G7 Product Catalog"*



*For diamond perfect performance*<sup>®</sup>

**accurate surgical & scientific instruments corporation**

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www accuratesurgical.com](http://www accuratesurgical.com)



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
*For diamond perfect performance*®