

ASSI[®] Plastic Surgery Instruments

ASSI UNDERMINER
ASSI.ATK50326

assi[®] 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS[®]

For diamond perfect performance[®]

accurate surgical & scientific instruments corporation
300 Shames Drive, Westbury, NY 11590
800.645.3569 516.333.2570 fax: 516.997.4948
west coast: 800.255.9378 • Info: assi@accuratesurgical.com
Orders: orders@accuratesurgical.com
www.accuratesurgical.com

UNDERMINERS



ASSI.ATK50326

- The NEW ASSI® Underminer is helpful in creating discontinuous tunnels of the lateral thigh and to immobilize the lateral thigh during the lower body lift procedure.
- The flat end is inserted along the deep fascia and several tunnels are created to enhance mobility of the lateral thigh.

ASSI® UNDERMINER

For blunt undermining in large procedures

ASSI.ATK50326

38cm, 13"; lift 3cm



ASSI.37362

LAROE UNDERMINING FORCEPS

ASSI.37362

24.5cm, 9³/₄", curved

ASSI.42362

30cm, 12¹/₂", curved

ASSI.52362

30cm, 12¹/₂", full curve

DEMARCATORS



ASSI.ATK39226

MARTEN FLAP FORCEPS

For accurate marking of facial flap

ASSI.ATK39226

14.5cm, 5³/₄", without ratchet lock

ASSI.ATK39326

14.5cm, 5³/₄", with ratchet lock



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

fax: 516.997.4948

800.645.3569 516.333.2570

west coast: 800.255.9378 • Info: assi@accuratesurgical.com

Orders: orders@accratesurgical.com

www accuratesurgical.com

DEMARCATORS

- The NEW ASSI® abdominal demarcator will help estimate the amount of tissue to be resected in body lift surgeries, especially abdominoplasty.
- Designed with an angle at the hinge, the ASSI abdominal demarcator accommodates the body's angle during surgery.
- Additionally, the longer and farther apart arms of the ASSI demarcator allow for thicker, larger flap resections.

ASSI.ATK391326



ASSI.ATK392326

ASSI ABDOMINAL DEMARCATORS

ASSI.ATK391326 small
160mm, 6 1/4"

ASSI.ATK392326 large
310mm, 12 1/2"



ASSI.ATK23826



ASSI.ATK389326

PITANGUY FLAP DEMARCATORS

ASSI ATK23826 small
15.5cm, 6 1/8"

ASSI ATK389326 large
26cm, 10 1/4"

Additional Demarcators on page 4

assi® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

fax: 516.997.4948

800.645.3569 516.333.2570


west coast: 800.255.9378 • Info: assi@accuratesurgical.com

Orders: orders@accuratesurgical.com

www accuratesurgical.com

DEMARCATORS

ASSI.ATK40926



ASSI.ATK64326

D'ASSUMPCAO FACE LIFT FORCEPS

For precise marking of dissected facial flaps for excision


ASSI.ATK40926
16cm, 6¹/₄", with ratchet

ASSI.ATK64326
16cm, 6¹/₄", without ratchet

ASSI.ATK41026
26cm, 10¹/₄", with ratchet

FACE LIFT FLAP FORCEPS

ASSI.ATK64526




ZINS FACE LIFT FLAP FORCEPS

ASSI.ATK64526
16cm, 6¹/₄", without ratchet

FLAP GRASPING FORCEPS


ASSI.ATK39026



PITANGUY FLAP GRASPING FORCEPS

ASSI.ATK39026
14cm, 5¹/₂"

ASSI.AG15326



LAHEY GRASPING FORCEPS

ASSI.AG15326
15cm, 6", 3x3 sharp teeth

assi® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

fax: 516.997.4948

800.645.3569 516.333.2570

west coast: 800.255.9378 • Info: assi@accuratesurgical.com

Orders: orders@accratesurgical.com

www accuratesurgical.com