

Surgical Instrument Color Coding Tapes



ASSI® Color Coding Tapes make it easy to create your own labeling system for identifying and sorting instruments by type, specialty, even ownership.

These durable tapes come in 10 high-visibility colors. They are easy to apply. They have been clinically tested and found safe for use in conventional, flash and gas sterilization units.

ASSI Color Coding Tapes are 1/4" wide and packaged on spools in 300" lengths. Enough tape to color code approximately 230 instruments.

ASSI.823	Green	
ASSI.825	Yellow	
ASSI.827	Blue	
ASSI.829	Red	
ASSI.831	Brown	
ASSI.833	White	
ASSI.835	Orange	
ASSI.837	Black	
ASSI.839	Purple	
ASSI.841	Glo-Orange	

Tip Covers – to prevent damage to instrument tips and edges

Now available vented or non-vented



- 11 Different sizes, each one color coded
- Autoclavable and steam-permeable
- Vented covers aid in the sterilizing and drying of condensation buildup on tips
- Stretchable — conforms to various instrument tip configurations
- Disposable
- Last for months
- Can be used on scissors, forceps, chisels, microsurgical and laparoscopic instruments, hemostats, curettes, rasps, knives, osteotomes, gouges, hooks, calipers, etc.
- Can be cut to any length
- Non-sterile

- ATC200 Mixed Assortment** (100 pieces, equal number of all sizes except ATC111-L)
- ATC200-V Mixed Assortment-Vented** (100 pieces, equal number of all sizes except ATC117 and ATC111-L)
- ATC400 Microsurgical Assortment** (100 pieces - 25 each: ATC100, 110, 111 and 112)
- ATC400-V Microsurgical Assortment-Vented** (100 pieces - 25 each: ATC100-V, 110-V, 111-V and 112-V)
- ATCBG Select Microsurgical Assortment** (100 pieces - 50 each: ATC100 and ATC110)
- ATCBG-V Select Microsurgical Assortment-Vented** (100 pieces - 50 each: ATC100-V and ATC110-V)
- Individual Package** (100 pieces, order by product number & color)
- Non-vented**
- Vented**

VENTED AND NON-VENTED SIZES AND COLORS											Tip covers shown actual sizes
NON-VENTED											
ATC000	ATC100	ATC110	ATC111	ATC112	ATC113	ATC114	ATC115	ATC116	ATC117	ATC111-L	
VENTED											
ATC000-V	ATC100-V	ATC110-V	ATC111-V	ATC112-V	ATC113-V	ATC114-V	ATC115-V	ATC116-V			
WHITE	BLUE	GREEN	RED	YELLOW	CLEAR	MAROON	ORANGE	BROWN	BLACK	BLUE4	RED
											Sold separately 100 piece pkgs.
											Not available in vented

Not all ASSI products shown in our literature or on our website are available for sale in Canada

accurate surgical & scientific instruments corporation

300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516-333-2570 • FAX: 1-516-997-4948 • WEST COAST SALES: SAN DIEGO, CA: 619-296-9945
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

1-800-645-3569 (NEW YORK OFFICE) **1-800-255-9378** (WEST COAST OFFICE)

Not all ASSI products shown in our literature or on our website are available for sale in Canada

ASSI[®] Single Use Products

TRUE BLUE[®] CONTRAST BACKGROUND MATERIAL

True Blue[®] microsurgical background material from ASSI[®] greatly enhances visibility in the surgical site. It can be trimmed to the size needed, slips in place easily and remains flat during the procedure. True Blue[®] is non-toxic, non-pyrogenic and not made with natural rubber latex. Available in either blue or yellow, with Smooth surface.

Order it by name – **True Blue[®]** will make a big difference.



ASSI.TB44S
Smooth Surface



ASSI.TBY44S
Smooth Surface

Catalog No.	Color	Texture	Size	Quantity
ASSI.TB44S	Blue	Smooth	1" x 2" Rectangle	Box of 12
ASSI.TBY44S	Yellow	Smooth	1" x 2" Rectangle	Box of 12
1 per sterile package – 12 per box				

Disposable Bipolar Pencils



6230BPD

18 Gauge Straight Bipolar Pencil (Green)

Serrated shaft provides safe grip with a fingertip platform for fast identification of the coagulation angle. Single piece molded construction prevents fluid short-out. Can be used in wet and dry fields. Packaged sterile 10 to a box. (☒sterile). 4³/₈" (11cm) long.



6232BPD

18 Gauge Curved Bipolar Pencil (Green)

Provides all the features and advantages of the straight pencil but with a 45° curved extended shaft for low manipulation in a sterile field. Packaged sterile 10 to a box. (☒sterile). 4³/₈" (11cm) long.



6234BPD

20 Gauge Straight Bipolar Pencil (Green)

Identical to the 18 gauge pencil. It may be used for coagulating choroidal and scleral bleeders. Packaged sterile 10 to a box. (☒sterile). 5" (12.5cm) long. Extended shaft is 1¹/₈" (29mm) long.



6236BPD

23 Gauge Tapered Tip Bipolar Pencil (Green)

Provides all the features of the 20 gauge straight pencil, but with a tapered point for easy coagulation entry of deep anterior segment bleeders. Packaged sterile 10 to a box. (☒sterile). 4³/₈" (11cm) long.



4230BPD

Jewelers Straight Tip Bipolar Forceps

Single use only. (☒sterile). Packaged 10 to a box. 4" (10.2cm) long.



2310BPD

Jewelers Curved Tip Bipolar Forceps

Single use only. (☒sterile). Packaged 10 to a box. 4" (10.2cm) long.



2320BPD

Coaptation Bipolar Forceps

Single use only. (☒sterile). Packaged 10 to a box. 4" (10.2cm) long



Disposable Sterile Bipolar Cord (One time use only)

12 foot (366cm) long Disposable Bipolar Cords are flexible, lightweight and will fit standard-type terminals. EtO sterilized.

227000C

*Free sample of 22-7000-C Disposable Cord is available upon request.

accurate surgical & scientific instruments corporation

1-800-645-3569 (NEW YORK OFFICE)

1-800-255-9378 (WEST COAST OFFICE)

Not all ASSi products shown in our literature or on our website are available for sale in Canada

NEW!

Breast Reduction Compression Device from ASSI®



The proven compression device that facilitates breast reduction surgery.

How to use:

1. Wrap the Breast Reduction Compression Device around the base of the breast and put the band into a loop.
2. Tighten the breast by contracting the band.
3. Release the Breast Reduction Compression Device by pulling back the small plastic tab at the head.

**STERILE – ITEM
SINGLE USE ONLY**



ASSI.BRC5
Standard size for an **approximate
500-800 gram reduction.**
5 pcs. per set

ASSI.BRC5XL
Extra large size for an **approximate
1300 gram reduction.**
5 pcs. per set

assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS CORPORATION
300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516.333.2570 • FAX: 1.516.997.4948 • WEST COAST SALES: SAN DIEGO, CA: 619.296.9945
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

Not all ASSI products shown in our literature or on our website are available for sale in Canada