







Tendon Retrievers, Grasping & Braiding Forceps and Tendon Strippers

TENDON RETRIEVERS


Rigid or Flexible Shaft, Large or Small Jaw. Stainless Steel Construction. Autoclavable.

<p>ASI3811</p>  <p>RIGID SHAFT</p>	<p>ASI3811 Rigid Shaft, Small Jaw Tip dia: 1mm 1 x 2 Teeth Jaw length: 8mm</p>  <p>ASI3813 Rigid Shaft, Large Jaw Tip dia: 1.5mm 1 x 2 Teeth Jaw length: 8mm</p> 
<p>ASI3812</p>  <p>FLEXIBLE SHAFT</p>	<p>ASI3812 Flexible Shaft, Small Jaw Tip dia: 1mm 1 x 2 Teeth Jaw length: 8mm</p> 
<p>ASSI.AG2668</p> 	<p>ASSI.AG2668 Carroll Tendon Retriever</p>



TENDON GRASPING FORCEPS

ASSI.AG1563




TENDON GRASPING FORCEPS WITH NEEDLE TEETH

Jaws feature two pairs of pin teeth that interdigitate upon closure. Limit to closure is adjustable by the thumb-screw attached to the instrument. Ratchet lock holds jaws firmly in place.

ASSI.AG1563
15cm (6")
Four pin teeth

ASSI.AG5339



TENDON GRASPING FORCEPS

Substantially reduce permanent tendon damage. Sharp teeth assure a solid grasp without tearing, even when applied to tendon ends.

ASSI.AG5339
15cm (6"), 2.5mm wide jaws
2x2 pin teeth

ASSI.AG1568




ISELIN TENDON GRASPING FORCEPS

Also known as Ombredanne Tendon Clamp, they feature "O" shaped jaws with 4x5 atraumatic teeth.

ASSI.AG1568
16cm (6-1/4")

ASSI.AG1569



ENGLERT TENDON GRASPING FORCEPS

Lightweight Englert forceps feature 4x5 atraumatic teeth.

ASSI.AG1569
14.5cm (5-3/4")

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS CORPORATION



300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516.333.2570 • FAX: 1.516.997.4948 • WEST COAST SALES: SAN DIEGO, CA: 619.296.9945
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com



Flip The Page
For More
INSTRUMENTS!



TENDON STRIPPERS

lightweight, autoclavable handle



 <p>3mm INNER DIA.</p>	ASSI.AG5331328 	3mm inner dia. ASSI.AG5331328 28.5cm (11") ASSI.AG1524335 35.0cm (14") ASSI.AG1525360 60.0cm (24")
---	---	---

 <p>4mm INNER DIA.</p>	ASSI.AG1524435 	4mm inner dia. ASSI.AG1524435 35.0cm (14") ASSI.AG1525460 60.0cm (24")
---	---	--



 <p>5mm INNER DIA.</p>	ASSI.AG5327528 	5mm inner dia. ASSI.AG5327528 28.5cm (11") ASSI.AG1524535 35.0cm (14") ASSI.AG5329550 50.0cm (20") ASSI.AG1525560 60.0cm (24")
---	---	--

 <p>6mm INNER DIA.</p>	ASSI.AG1524635 	6mm inner dia. ASSI.AG1524635 35.0cm (14") ASSI.AG1525660 60.0cm (24")
---	---	--


ADDITIONAL TENDON STRIPPERS CAN BE FOUND ON NEXT PAGE →

TENDON BRAIDING FORCEPS (TENDON WEAVING FORCEPS)


Provide valuable assistance during tendon braiding. The conically shaped fine point on one jaw separates tendon without causing long splits. The jaws can grasp and hold a small tendon stump securely.





ASSI.AG5343
15cm (6") **adult size**, curved
Delicate sharp point
Single pin tooth



ASSI.AG15742
12cm (4-3/4") **child size**, curved
Delicate sharp point
Single pin tooth




ASSI.AG5341
15cm (6") **adult size**, straight
Delicate sharp point
Single pin tooth



ASSI.AG15741
12cm (4-3/4") **child size**, straight
Delicate sharp point
Single pin tooth

BRAND TENDON TUNNELING/PASSING FORCEPS






ASSI.AG79026
25cm (9-7/8")

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS CORPORATION

300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516.333.2570 • FAX: 1.516.997.4948 • WEST COAST SALES: SAN DIEGO, CA: 619.296.9945
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com


HALSTED-S TENDON GRASPING FORCEPS

These forceps simplify the task of maneuvering a tendon around a phalanx.

<p>ASSI.AG15711</p>  <p>A pair of stainless steel tendon grasping forceps with a light curve. The handles are oval-shaped and the jaws are long and thin with a slight upward curve. A circular inset shows a magnified view of the serrated, curved jaws.</p>	<p>ASSI.AG15711 12.5cm (5") Light curve</p>
<p>ASSI.AG15712</p>  <p>A pair of stainless steel tendon grasping forceps with a medium curve. The handles are oval-shaped and the jaws are long and thin with a moderate upward curve. A circular inset shows a magnified view of the serrated, curved jaws.</p>	<p>ASSI.AG15712 12.5cm (5") Medium curve</p>
<p>ASSI.AG15713</p>  <p>A pair of stainless steel tendon grasping forceps with a strong curve. The handles are oval-shaped and the jaws are long and thin with a pronounced upward curve. A circular inset shows a magnified view of the serrated, curved jaws.</p>	<p>ASSI.AG15713 12.5cm (5") Strong curve</p>

DIETRICH TENDON FIXATION NEEDLE

Temporarily immobilizes the tendon, restricting muscular action during the braiding process.

<p>ASSI.AG5345</p>  <p>A long, thin, tapered stainless steel needle with a rounded handle. The needle is designed for tendon fixation during braiding.</p>	<p>ASSI.AG5345 4.5cm (1-3/4")</p>
--	--










ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS CORPORATION

300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516.333.2570 • FAX: 1.516.997.4948 • WEST COAST SALES: SAN DIEGO, CA: 619.296.9945
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

TENDON STRIPPERS

BUNNELL TENDON STRIPPER









13.5cm (5 1/2")

<p>ASSI.AG152725</p> 	 2.5mm dia.	<p>ASSI.AG152725 2.5mm dia.</p>
<p>ASSI.AG152730</p> 	 3.0mm dia.	<p>ASSI.AG152730 3.0mm dia.</p>
<p>ASSI.AG152740</p> 	 4.0mm dia.	<p>ASSI.AG152740 4.0mm dia.</p>
<p>ASSI.AG152750</p> 	 5.0mm dia.	<p>ASSI.AG152750 5.0mm dia.</p>
	 6.5mm dia.	<p>ASSI.AG152765 6.5mm dia.</p>

FLEXIBLE SHAFT TENDON STRIPPER

22cm (8 3/4")



<p>ASSI.AG152840</p> 	 4.0mm dia.	<p>ASSI.AG152840 4.0mm dia.</p>
<p>ASSI.AG152855</p> 	 4.5mm dia.	<p>ASSI.AG152845 4.5mm dia.</p>
	 5.5mm dia.	<p>ASSI.AG152855 5.5mm dia.</p>
<p>ASSI.AG152860</p> 	 6.0mm dia.	<p>ASSI.AG152860 6.0mm dia.</p>
	 7.0mm dia.	<p>ASSI.AG152870 7.0mm dia.</p>

accurate surgical & scientific instruments corporation

300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516.333.2570 • FAX: 1.516.997.4948 • WEST COAST SALES: SAN DIEGO, CA: 619.296.9945

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

1.800.645.3569 (NEW YORK OFFICE) **1.800.255.9378** (WEST COAST OFFICE)