

# Unique Engler Retractor For Facial Procedures


Provides Maximum Exposure  
and Vision for Small Incisions.



Submental Incision, showing reach of retractor



Submental Incision, retractor in place



- Ideal for Liposuction of the Chin and Neck, Facelift, Browlift, Facial Implants, Suspension Devices and more
- Narrow width of blade makes it ideal for use with the smallest incisions
- Long length of blade delivers light and retraction much farther from the incision than standard retractors
- Compatible with all standard, existing light sources and fiber optic cables

Designed by Alan Engler, M.D.

Patented


ASSI.ABR54326  
65mmx12mm wide blade  
with Fiber Optic and Adapter

ASSI.ABR54226  
65mmx12mm wide blade  
with Suction  
w/o Fiber Optic

ASSI.ABR54726  
65mmx12mm wide blade, **insulated**,  
with Fiber Optic and Adapter

ASSI.ABR54626  
65mmx12mm wide blade

ASSI.ABR54026  
65mmx12mm **insulated blade with  
teeth**, Fiber Optic and Adapter



ASSI.ABR54426  
65mmx12mm wide blade **with teeth**,  
Fiber Optic and Adapter

ASSI.ABR54526  
65mmx12mm wide blade **with teeth**

**assi**®   
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®

For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

[www.accuratesurgical.com](http://www.accuratesurgical.com)