


# ASSI® Endoscopes for Plastic Surgery

- Ideal for procedures such as Face Lifts, Forehead Lifts, Breast Augmentation and Abdominoplasty
- High quality German optics
- Superior visualization in smaller pockets (4 & 5mm scope) as well as in larger pockets (10mm scope)
- The 30° lens angle allows for rotation of the scope, increasing the field of vision without repositioning the scope



**ASSI.AB382126**  
4mm dia., 0° wide angle, 17.5cm working length

**ASSI.AB382226**  
5mm dia., 0° wide angle, 17.5cm working length

**ASSI.AB383126**  
4mm dia., 30° wide angle, 17.5cm working length


**ASSI.AB383226**  
5mm dia., 30° wide angle, 17.5cm working length



**ASSI.AB384126**  
10mm dia., 30° wide angle, 17.5cm working length

**ASSI.AB384226**  
10mm dia., 0° wide angle, 17.5cm working length

**ASSI.AB384326**  
10mm dia., 30° wide angle, with Storz adapter, 17.5cm working length



**AS3293C**  
See-Thru Medium Scope Case,  
15.5 x 2.75 x 1.38 inch,  
Holds 1 Scope 1-10mm Dia, Up to 15.5 inch long

**AS3323C**  
See-Thru Large Scope Case,  
18 x 2.75 x 1.38 inch,  
Holds 1 Scope 1-10mm Dia, Up to 18 inch long

**assi®**   
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
For diamond perfect performance®

accurate surgical & scientific instruments corporation  
300 Shames Drive, Westbury, NY 11590  
800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378  
Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)  
[www.accuratesurgical.com](http://www.accuratesurgical.com)