

ASSI Gutow Self Retaining Retractor

For Use In Hand Surgery.



ASSI.ATK17826
Gutow Self Retaining Retractor

- The Gutow Self Retaining Retractor is ideally designed for Minimal Incision Dorsal Scaphoid Fixation. The blunt tips are precisely configured for intra-articular retraction.
- Additionally, it is useful for common hand surgery procedures such as Open Ganglion Excisions where the blunt tips can retract the surrounding tendons and vessels without injury.
- The Gutow Self Retaining Retractor is also uniquely useful for Trigger Finger Release. The flat blunt retractors can delicately retract the neurovascular bundles without damage.

Designed by:
Andrew Gutow, M.D. FAAOS, Palo Alto, CA

assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590 • 800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

Not all ASSI products shown in our literature or on our website are available for sale in Canada