





ASSI® Levine Breast Retractor & Implant Protector

Features:

- Designed to facilitate closure of the tissues over a prosthesis
- This Retractor simplifies exposure and eliminates accidental implant damage when suturing acellular matrix in tissue expander or immediate breast reconstruction.
- Prevents accidental damage to saline filled implants
- The Levine Breast Retractor & Implant Protector is useful in all breast prosthetic procedures including augmentations as well as reconstructions.



ASSI.AB21426
Levine Breast Retractor
& Implant Protector
Double Ended

Designed by:
Richard Levine, MD
San Antonio, Texas

assi® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

Not all ASSI products shown in our literature or on our website are available for sale in Canada