

ASSI® Micro Monopolar Forceps

*Insulated, hand-activated
non-stick coating, straight, smooth jaws.*

ASSI® Micro Monopolar Forceps are specifically constructed to reach the entire length of the pocket, and as such are ideal for inframammary and initial, upper axillary pocket dissections. The forceps precisely combines cut and coagulation current for best possible dissection, while averting unnecessary current spread, tissue damage and bleeding. Limited reuse.

401500MMP

15cm, 6 in., 0.7mm tip dia.

bayonet

401900MMP

19cm, 7½ in., 0.7mm tip dia.

bayonet

401240MMP

24cm, 9½ in., 0.7mm tip dia.

bayonet

402400MMP


24cm, 9½ in., 0.7mm tip dia.

straight

401270MMP

26cm, 10¼ in., 0.7mm tip dia.

straight



assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com
www.accuratesurgical.com