


Unique Engler Retractor For Breast Augmentation


Provides Maximum Exposure
and Vision for Small Incisions.



Periareolar Incision, showing reach of retractor



Periareolar Incision, retractor in place



- Ideal for the Periareolar, Inframammary and Trans-axillary incisions
- Narrow width of blade makes it ideal for use with the smallest incisions
- Long length of blade delivers light and retraction much farther from the incision than standard retractors
- Blade provides exposure for both dissection and hemostasis
- Compatible with all standard, existing light sources and fiber optic cables

Designed by Alan Engler, M.D.
Patented


ASSI.ABR54326
65mmx12mm wide blade
with Fiber Optic and Adapter

ASSI.ABR54226
65mmx12mm wide blade
with Suction
w/o Fiber Optic

ASSI.ABR54726
65mmx12mm wide blade,
insulated,
with Fiber Optic and Adapter

ASSI.ABR54626
65mmx12mm wide blade
w/o Fiber Optic and Suction

ASSI.ABR54026
65mmx12mm wide **insulated** blade
with teeth, Fiber Optic and Adapter



ASSI.ABR54426
65mmx12mm wide blade **with teeth**,
Fiber Optic and Adapter

ASSI.ABR54526
65mmx12mm wide blade **with teeth**
w/o Fiber Optic

assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation
300 Shames Drive, Westbury, NY 11590
800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com
www.accuratesurgical.com