

# CLEANING BRUSHES

## FOR LARGE SURGICAL INSTRUMENTS


Proper cleaning can substantially increase the life of your surgical instruments. Soaking in cleaning solution alone doesn't always get the job done. That's because an instrument's small serrations and crevices can trap particles that require hand scrubbing.

Our specialty brushes are the answer.

They're perfect for precleaning serrated surfaces and box locks before sterilization in the Autoclave. The more delicate nylon bristles can be safely used on any instrument surface. The more rugged stainless steel bristles are recommended for brushing raps, files and burrs.

Why not take every precaution? Make these cleaning brushes a part of your daily maintenance routine.

**AMS207**  
Stainless steel  
bristle cleaning brush



**AMS205**  
Nylon bristle  
cleaning brush

CAT. NO.	DESCRIPTION
AMS205	Nylon bristle cleaning brush
AMS207	Stainless steel bristle cleaning brush
Set contains 3 of the same type of brush	



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
*For diamond perfect performance®*

**accurate surgical & scientific instruments corporation**

300 SHAMES DRIVE, WESTBURY, NY 11590, TEL: 516.333.2570 • FAX: 1.516.997.4948 • WEST COAST SALES: 619.296.9945

**1.800.645.3569** (NEW YORK) **1.800.255.9378** (WEST COAST)

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)