


ASSI[®] Perma Facial Implant[™] Tunneler Passers and Forceps




Perma Implant
Tunneler
Passer Small
ASSI.80320



Perma Implant
Tunneler
Passer Large
ASSI.80326



Perma Implant
Forceps
Straight
ASSI.78026



Perma Implant
Forceps
Angled
ASSI.77926

assi 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS[®]
For diamond perfect performance[®]

accurate surgical & scientific instruments corporation
516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
800.645.3569 • Info: assi@accuratesurgical.com
Orders: orders@accuratesurgical.com
www.accuratesurgical.com



ASSI Perma Implant Tunneler Passer Small with Perma Facial Implant



ASSI Perma Implant Tunneler Passer Small, Perma Facial Implant and Perma Implant Forceps Straight



ASSI Perma Implant Tunneler Passer Large with Perma Facial Implant




ASSI Perma Implant Tunneler Passer Large, Perma Facial Implant and Perma Implant Forceps Angled



Perma Facial Implant is a trademark of SurgiSil L.L.P.

ASSI® Perma Facial Implant™ Tunnelers Passers and Forceps



The ASSI® Perma Implant Tunnelers and Perma Implant Forceps were specifically designed to be used with the Perma Facial Implants. The Large and Small ASSI Implant Tunneler Passers are ideally suited for passage through the subcutaneous tissue or lip, and securely grasp the Perma Facial Implant. The Large ASSI Implant Tunneler Passer is better for larger implants, and the smaller ASSI Implant Tunneler Passer is better for the smaller Perma Facial Implant. The unique design of the jaws on both the ASSI Implant Tunneler Passers and ASSI Perma Implant Forceps facilitates proper grasping of the Perma Facial Implant without causing trauma to the implant.




ASSI.80320
Perma Implant Tunneler
Passer Small
15cm, special oval jaws TCC



ASSI.80326
Perma Implant Tunneler
Passer Large
19.5cm, special oval jaws TCC



ASSI.78026
Perma Implant Forceps
Straight
12cm, special oval jaws TCC



ASSI.77926
Perma Implant Forceps
Angled
12cm, special oval jaws TCC