




ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®

# Lalonde® Percutaneous Bone Clamp with K-Wire Guide

## “Finally, A Percutaneous Bone Clamp For Finger Fractures . . . That Works!”

**Designed For Closed Reduction And Percutaneous K-Wire Fixation Of Phalanx Fractures.**

- This finger fitting bone clamp has sharp tips that reduce and hold fractures like a towel clip.
- The K-Wire guide on top of the instrument aims the K-Wire right where you want it — aims the K-Wire through the fracture to just above the reducing tips.




ASSI.AG776626

Lalonde® Percutaneous Bone Clamp With K-Wire Guide  
Max. K-Wire .045", 12.5cm

Designed by D.H. Lalonde, B.Sc., M.Sc., M.D., FRCSC  
Saint John, N.B., Canada

Patented



ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®

For diamond perfect performance®

accurate surgical & scientific instruments corporation

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)

[www.accuratesurgical.com](http://www.accuratesurgical.com)

Not all ASSI products shown in our literature or on our website are available for sale in Canada