

# Di Geronimo Breast Clamp

ASSI-AB80326



## Fill Tube Inserter Clamp

- Use to reinsert the fill tube into the breast implant valve.
- Saves time.
- Special jaw configuration securely holds fill tube end, making reinsertion easy.
- Use if the tube slips out during filling procedure.
- Use if more or less saline solution is required.
- Decreases risk to breast implant integrity.

**assi**®   
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
*For diamond perfect performance®*

accurate surgical & scientific instruments corporation  
300 Shames Drive, Westbury, NY 11590  
800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378  
Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com)  
[www.accuratesurgical.com](http://www.accuratesurgical.com)