

Newman Collagen Plug Inserter Forceps

FOR PRECISE CONTROL OF COLLAGEN PLUGS

The Newman Collagen Plug Inserter Forceps features a specially designed channel at the tip so it is easy to grasp the plug out of the kit and a special locking mechanism that is easy to open and close with the thumb or index finger. (See figure A)

Unlike jewelers forceps, which have only two-point fixation, this unique channel design minimizes the risk of the collagen plug slipping out or being crushed. (See figure B)

These features make the Newman Collagen Plug Inserter Forceps a more efficient and cost effective instrument.

ASSI.260529
11cm, with lock, right

ASSI.261529
11cm, with lock, left




FIGURE A




FIGURE B

Designed by Bruce Newman, M.D.
Phoenix, AZ
Patent Pending

assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation
300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com