


Campbell Lip Augmentation Awl

For creating a submucosal tunnel;
insertion and placement of graft.



ASSI.38362

Curved, with suture hole


54mm overall length




*Close-up detail of
suture hole at tip*



*Making horizontal
incisions with ASSI®
Campbell Lip
Augmentation Awl*



*Tunneling with ASSI
Campbell Lip
Augmentation Awl*



*Pulling graft
material through
tunnel with the
ASSI Campbell Lip
Augmentation Awl*

*Enlarged to show
detail*

assi® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378 www accuratesurgical.com