



Lipshultz Epididymovasostomy Microdissection Scissors

“These unique scissors enable the careful dissection not only of the tunica of the epididymis for epididymal tubular exposure, but also allow elliptical cutting of the tubule itself to ensure a perfect circle for an epididymovasostomy reanastomosis.”



Larry I. Lipshultz, M.D.,
Professor of Urology,
Baylor College of Medicine



ASSI.SDC15RVL

- Decreases the chance of inadvertently opening an epididymal tube during dissection.
- Well-balanced, resting easily in the hand.
- Unusually smooth ventrally.

ASSI.SDC15RVL

Dissecting Scissors, Lipshultz pattern, 15cm, round handle, 8mm diameter, **curved, slender, tapered**, 7mm long blades

ASSI.SAS15RVL

Dissecting Scissors, Lipshultz pattern, 15cm, round handle, 8mm diameter, **straight, slender, tapered**, 7mm long blades



accurate surgical & scientific instruments corporation
300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com