




# NO-SCALPEL VASECTOMY INSTRUMENTS

## NO-SCALPEL VASECTOMY FIXATOR RING CLAMP FORCEPS




**ASSI.VAS47526**  
14cm, 5.5 in.; with blunt tips, **standard ring**, Gold handles




**ASSI.VAS46326**  
14cm, 5.5 in.; with blunt tips, **small ring**, Gold handles


## NO-SCALPEL VASECTOMY FORCEPS



**ASSI.VAS47626**  
12.5cm, 5 in.; with sharp pointed smooth curved jaw, Gold handles




**ASSI.VAS47726**  
12.5cm, 5 in.; **very delicate curved**, smooth pointed jaw, Gold handles



## NO-SCALPEL VASECTOMY INSTRUMENT SETS




**ASSI.VAS 95126**  
*Consists of:*  
1 each ASSI.VAS47626 Forceps  
1 each ASSI.VAS47526 Fixator Ring Clamp, **standard ring**




**ASSI.VAS 93952**  
*Consists of:*  
1 each ASSI.VAS47626 Forceps  
1 each ASSI.VAS46326 Fixator Ring Clamp, **small ring**

**1** Fixation of the vas under the scrotal skin using the extracutaneous vas fixation clamp, ASSI® No-Scalpel Vasectomy Fixator Ring Clamp.



**2** Puncture of the scrotal skin using the ASSI® No-Scalpel Vasectomy Forceps with sharp pointed smooth curved jaw.



**3** Dissecting blood vessels and sheath of vas prior to division and ligation or cauterization of the lumen.

