



ASSI Haynes Breast Retractor




- The proximal teeth help to hold the tissue on tension while scoring with the cautery



- Permits the surgeon to easily modify an implant pocket – Especially scoring the capsule to adjust the shape or treat contracture



- The slotted blade allows the surgeon to use the cautery to cut up close to the base of the Retractor




Designed By:
Daniel Haynes M.D. Johnson City, TN

ASSI.ABR83826

Haynes Breast Retractor
90 x 32mm slotted blade,
with teeth (Shown)

ASSI.ABR84226

Haynes Breast Retractor
140 x 32mm slotted blade,
with teeth



ASSI.ABR84026

Haynes Breast Retractor
90 x 32mm slotted blade,
with teeth, w/fiber optic

ASSI.ABR84426

Haynes Breast Retractor
140 x 32mm slotted blade,
with teeth, w/fiber optic (Shown)



assi® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com

www.accuratesurgical.com

Not all ASSI products shown in our literature or on our website are available for sale in Canada