


TIP COVERS

Designed to prevent tip & edge damage to instruments



Now Available
Non-vented or Vented
100 pieces per pack

- 11 Different sizes, each one color coded
- Autoclavable and steam-permeable
- Vented covers aid in the sterilizing and drying of condensation buildup on tips
- Stretchable — conforms to various instrument tip configurations
- Disposable
- Last for months
- Can be used on scissors, forceps, chisels, microsurgical and laparoscopic instruments, hemostats, curettes, rasps, knives, osteotomes, gouges, hooks, calipers, etc.
- Can be cut to any length
- Non-sterile

ATC200 Mixed Assortment
(100 pieces, equal number of all sizes except ATC111-L)

ATC200-V Mixed Assortment-Vented
(100 pieces, equal number of all sizes except ATC117 and ATC111-L)

ATC400 Microsurgical Assortment
(100 pieces - 25 each: ATC100, 110, 111 and 112)

ATC400-V Microsurgical Assortment-Vented
(100 pieces - 25 each: ATC100-V, 110-V, 111-V and 112-V)

ATCBG Select Microsurgical Assortment
(100 pieces - 50 each: ATC100 and ATC110)

ATCBG-V Select Microsurgical Assortment-Vented
(100 pieces - 50 each: ATC100-V and ATC110-V)

Individual Package (100 pieces, order by product number & color)
Non-vented
Vented

VENTED AND NON-VENTED SIZES AND COLORS

Tip covers shown actual sizes

NON-VENTED

ATC000 ATC100 ATC110 ATC111 ATC112 ATC113

ATC114

ATC115

ATC116

ATC117

ATC111-L

VENTED

ATC000-V ATC100-V ATC110-V ATC111-V ATC112-V

ATC113-V

ATC114-V

ATC115-V

ATC116-V

Sold separately
100 piece pkgs.



WHITE



BLUE



GREEN



RED



YELLOW



CLEAR



MAROON



ORANGE



BROWN



BLACK



BLUE4



RED

Not available in vented

accurate surgical & scientific instruments corporation

1.800.645.3569 (NEW YORK OFFICE)

1.800.255.9378 (WEST COAST OFFICE)

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com