


assi®

AUFRICHT NASAL RETRACTOR

INSULATED WITH PERFORATED SUCTION TUBE

as designed by Dr. Theodore Diktaban

- For use with **VAPOR** (*Vaporized Assisted Profile Option Reduction*) **Nose Technique** developed by Dr. Theodore Diktaban.
- Enhanced visualization for closed rhinoplasty using **VAPOR Nose Technique**.
- Enhanced safety with the insulated working end which decreases heat conduction to the overlying skin during the **VAPOR Nose Technique**.
- Evacuates plume and maintains visualization during the **VAPOR Nose Technique**.
- Can be used in open rhinoplasty procedures when the electrocautery is used.



ASSI.AG57326

DESIGNED BY:
Dr. Theodore Diktaban
New York, NY

assi® 
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation
300 Shames Drive, Westbury, NY 11590
800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378
Info: assi@accuratesurgical.com
Orders: orders@accuratesurgical.com • www.accuratesurgical.com
Note: Not all ASSI products shown in our literature or on our website are available for sale in Canada.

©2020 ASSI®
Spring-4M-2MFT.2MFD-6
F-168