


# ASSI®



Designed By:  
David A. Hidalgo, M.D.  
Clinical Professor of Surgery  
Cornell University Medical College  
New York, NY

## HIDALGO™ SUBMENTAL FACE LIFT RETRACTOR

### Ideal Exposure For:

- Caudal aspect of platysmal plication.
- Low transverse platysmal division.
- Anterior cervical defatting and hemostasis.
- Medial SMAS dissection and fixation.

### Retractor Features:

- Curved side wings expand visualization
- Fiberoptic endpoint elongates the optical cavity
- Polished Surfaces enhance illumination

ASSI.ABR90626

55 x 12mm wide blade w/ Wings, Fiber Optics & Universal Adapter

ASSI.ABR90426

55 x 12mm wide blade w/ Wings, Without Fiber Optics & Universal Adapter



ASSI.ABR90626

**assi®**   
ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®  
*For diamond perfect performance®*

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590 • 800.645.3569

ph: 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: [assi@accuratesurgical.com](mailto:assi@accuratesurgical.com) • Orders: [orders@accuratesurgical.com](mailto:orders@accuratesurgical.com) • [www.accuratesurgical.com](http://www.accuratesurgical.com)

NOTE: Not all ASSI products shown in our literature or on our website are available for sale in Canada