




ASSI[®] BREAST

RETRACTOR WITH WEIGHT

“The Hands-Free Breast Deaver”

- Frees up the Assistant's hand
- Eliminates the cumbersome need to reach across the OR table
- The curved, channeled blade tip grasps the tissue edge to hold the retractor in position
- The weighted end becomes an ergonomic handle when turned upwards



Designed By:

Kenneth R. Francis, MD, FACS
Clinical Assistant Professor of Plastic Surgery
The New York Presbyterian Hospital
Cornell Weill Medical College
New York, NY

ASSI.AG73126

assi[®] 
DFFXUDWH#XUJ IEFDO# #F#HQLV IEF#QVWXP HQWV[®]
For diamond perfect performance[®]

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590 • 800.645.3569

ph: 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com • www.accuratesurgical.com

NOTE: Not all ASSI products shown in our literature or on our website are available for sale in Canada