

ASSI®'s Lipshultz Urology Microsurgical Set

For Vasovasostomy and Vasoepididymostomy

Recommended by Larry I. Lipshultz, M.D.,
Professor of Urology,
Baylor College of Medicine

ASSI.LVAS1 SET Lipshultz Urology Microsurgical Set

Set features Lipshultz Epididymovasostomy Dissecting Scissors, including instruments tailored specifically for Vasovasostomy and Vasoepididymostomy surgery.



Set consists of the following: (no substitutions)

QUANTITY	CATALOG NUMBER	PRODUCT DESCRIPTION
Four	ASSI.4318	Castroviejo Forceps, 12mm
One	ASSI.B158	Round Handle Needle Holder
One	ASSI.SDC15RVL	"Lipshultz" Epididymovasostomy Dissecting Scissors
One	103000BP	Microtip Bipolar Forceps
Two	ASIM5511	Jewelers Forceps, 0.3mm Tip
Two	ASSI.JF4	Jewelers Forceps, 0.4mm Blunt Tip
One	ASSI.NHF2	Vas Holding Forceps, 2mm
One	ASSI.NHF2.5	Vas Holding Forceps, 2.5mm
One	ASSI.NHF3	Vas Holding Forceps, 3mm
One	ASSI.NHF3.5	Vas Holding Forceps, 3.5mm
One	ASSI.BHS11	Blade Holder, 11cm
One Box	ASSI.CBS35	Straight Blades, Sterile (10 blades/box)
One	AS1174C	Instrument Case with Silicone Mat

assi® 

ACCURATE SURGICAL & SCIENTIFIC INSTRUMENTS®
For diamond perfect performance®

accurate surgical & scientific instruments corporation

300 Shames Drive, Westbury, NY 11590

800.645.3569 516.333.2570 fax: 516.997.4948 west coast: 800.255.9378

Info: assi@accuratesurgical.com • Orders: orders@accuratesurgical.com

www.accuratesurgical.com

Not all ASSI products shown in our literature or on our website are available for sale in Canada